

## DW-VET13 Color Ultrasonic Diagnostic Apparatus



## Technical Specification

|           |                                                                                                                                                                                                                       |
|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>1:</b> | <b>Summary of main specifications and system of cart type color Doppler ultrasound</b>                                                                                                                                |
| 1.1       | Portable type all digital color Doppler ultrasonic mainframe                                                                                                                                                          |
| 1.2       | Ultrasonic host operating system: Windows operating system                                                                                                                                                            |
| 1.3       | Applications:<br>Abdomen, obstetrics, gynecology, heart, urinary system, small organs, superficial, blood vessels, pediatrics, newborns, musculoskeletal                                                              |
| 1.4       | Probes:<br>Convex probe, Tran-vaginal probe, Linear probe, Micro-convex probe, Cardiac probe, 4D Volume probe                                                                                                         |
| 1.5       | Applications and report:<br>Abdominal,OB,GYN,Cardiac,Urinary,Small Parts,Superficial, Vascular, Pediatrics, Advanced measurement software packages, report software packages, case management software packages, etc. |
| ☆1.6      | carotid artery intima measurement thickness(IMT)                                                                                                                                                                      |
| ☆1.7      | Automatic spectral envelope measurement                                                                                                                                                                               |
| 1.8       | Full digital transmission and reception of beam synthesizer                                                                                                                                                           |
| 1.9       | Color Doppler imaging (C)                                                                                                                                                                                             |
| 1.10      | Pulse Doppler Imaging (PW)                                                                                                                                                                                            |
| ☆1.12     | Coherent Contrast imaging (CCI)                                                                                                                                                                                       |
| 1.12      | Continuous wave Doppler imaging (CW)                                                                                                                                                                                  |
| ☆1.13     | B/C/D Real-time three synchronous imaging                                                                                                                                                                             |
| ☆1.14     | Power Doppler imaging (PDI)                                                                                                                                                                                           |



DAWEI MEDICAL (JIANGSU) CORP., LTD.

|                        |                                                                                                                                                                                                                             |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ☆1.15                  | Direct power Doppler imaging (DPDI)                                                                                                                                                                                         |
| 1.16                   | M mode imaging                                                                                                                                                                                                              |
| ☆1.17                  | Anatomic M mode imaging                                                                                                                                                                                                     |
| ☆1.18                  | Color Doppler M mode imaging                                                                                                                                                                                                |
| ☆1.19                  | Elastography                                                                                                                                                                                                                |
| ☆1.20                  | Tissue Doppler imaging (TDI)                                                                                                                                                                                                |
| ☆1.21                  | Strain rate imaging (SRI)                                                                                                                                                                                                   |
| 1.22                   | Tissue harmonic imaging (THI)                                                                                                                                                                                               |
| 1.23                   | Fusion harmonic imaging (FHI)                                                                                                                                                                                               |
| 1.24                   | Speckle Reduce imaging (SRI)                                                                                                                                                                                                |
| ☆1.25                  | Panoramic imaging                                                                                                                                                                                                           |
| ☆1.26                  | Deflection imaging                                                                                                                                                                                                          |
| ☆1.27                  | Trapezoidal imaging                                                                                                                                                                                                         |
| 1.28                   | Adaptive velocity optimization                                                                                                                                                                                              |
| ☆1.29                  | Free hand 3D                                                                                                                                                                                                                |
| 1.30                   | Real time 3D imaging (3D/4D)                                                                                                                                                                                                |
| 1.31                   | DICOM3.0                                                                                                                                                                                                                    |
| 1.32                   | Monitor: 15 inch, high definition ultrasonic display                                                                                                                                                                        |
| 1.33                   | Physical clipboard: save the image on the left side of the screen, which can be directly saved or deleted.                                                                                                                  |
| 1.34                   | The system has the function of on-the-spot upgrade                                                                                                                                                                          |
| 1.35                   | Presupposition: for different inspection of the viscera, preset the inspection conditions for the best image, reduce the adjustment of the operation, and the commonly used external adjustment and combination regulation. |
| 1.36                   | Probe interface: 2                                                                                                                                                                                                          |
| 1.37                   | Language: CN/EN/VI/DE/FR/ES/RU/AR/PT/ID                                                                                                                                                                                     |
| 1.38                   | Depth: ≥360mm;                                                                                                                                                                                                              |
| 1.39                   | Extended imaging                                                                                                                                                                                                            |
| <b>2:</b>              | <b>Probes</b>                                                                                                                                                                                                               |
| 2.1<br>Convex<br>probe | Fundamental Frequency: 2.0MHz/2.3MHz/2.5MHz/3.0MHz/3.5MHz/4.0MHz/4.6MHz/5.0MHz/5.4MHz,<br>Harmonic Frequency:<br>4.0MHz/4.6MHz/5.0MHz,                                                                                      |
| 2.2<br>Linear<br>probe | Fundamental Frequency:<br>4.0MHz/4.6MHz/5.0MHz/6.0MHz/7.0MHz/8.0MHz/9.2MHz/10.0MHz/12.0MHz/13.3MHz,<br>Harmonic Frequency:<br>8.0MHz/9.2MHz/10.0MHz,                                                                        |
| 2.3<br>Trans-vagin     | Fundamental Frequency:<br>3.0MHz/3.5MHz/4.0MHz/5.0MHz/5.4MHz/6.0MHz/7.0MHz/8.0MHz/10.0MHz,                                                                                                                                  |



DAWEI MEDICAL (JIANGSU) CORP., LTD.

|                               |                                                                                                                                                   |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| al probe                      | Harmonic Frequency:<br>6.0MHz/7.0MHz/8.0MHz,                                                                                                      |
| 2.4<br>Micro-conv<br>ex probe | Fundamental Frequency:<br>3.0MHz/3.5MHz/4.0MHz/5.0MHz/5.4MHz/6.0MHz/7.0MHz/8.0MHz,<br>Harmonic Frequency:<br>6.0MHz/7.0MHz/8.0MHz,                |
| 2.5<br>Cardiac<br>probe       | Fundamental Frequency: 1.7MHz/1.9MHz/2.1MHz/2.5MHz/3.0MHz/3.4MHz/3.8MHz/4.2MHz/5.0MHz,<br>Harmonic Frequency:<br>3.4MHz/3.8MHz/4.2MHz,            |
| 2.6<br>4D Volume<br>probe     | Fundamental Frequency:<br>2.0MHz/2.5MHz/3.0MHz/3.3MHz/3.7MHz/4.0MHz/5.0MHz/6.0MHz,<br>Harmonic Frequency:<br>4.0MHz/5.0MHz/6.0MHz,                |
| <b>3:</b>                     | <b>2D imaging mode</b>                                                                                                                            |
| 3.1                           | Gain: 0—100, Step 2 adjustable                                                                                                                    |
| 3.2                           | TGC: 8 segment adjustable                                                                                                                         |
| 3.3                           | Maximum focus point: $\geq 7$ , which can be moved throughout the whole process.                                                                  |
| 3.4                           | Speckle reduction: 0-5, 5 level                                                                                                                   |
| 3.5                           | Space Synthesis: 0-2, 2 level (Liner probe: 3 level, cardiac probe:0)                                                                             |
| 3.6                           | Dynamic: 30-180, 35 level, step 5 adjustable                                                                                                      |
| 3.7                           | Line density: low、middle、high, 3 level                                                                                                            |
| 3.8                           | Frame correlation: 0-4, 4 level                                                                                                                   |
| 3.9                           | Noise reduction: 0-5, 5 level                                                                                                                     |
| 3.10                          | Edge Enhancement: 0-5, 5 level                                                                                                                    |
| 3.11                          | Sound power: 2-10, 9 level                                                                                                                        |
| 3.12                          | Grey scale: 0-67, 67 level                                                                                                                        |
| 3.13                          | False color: 0-67, 67 level                                                                                                                       |
| 3.14                          | Image style: Soft-Comparison, 2 level                                                                                                             |
|                               | The screen has real-time display of voice power, probe frequency, dynamic range, pseudo color, gray scale and other 11 parameters can be adjusted |
| <b>4:</b>                     | <b>Color Doppler imaging mode</b>                                                                                                                 |
| 4.1                           | Blood gain: 0-100, Step 2                                                                                                                         |
| 4.2                           | Parameter display: Velocity、Variance                                                                                                              |
| 4.3                           | B-Restrain (B/W restrain) : 0-7, 7 level                                                                                                          |
| 4.4                           | Speed Through: 0-8, 8 level                                                                                                                       |
| 4.5                           | Sampling number: 6-24, 7 level                                                                                                                    |
| 4.6                           | Blood flow preferred: 0-8, 8 level                                                                                                                |

|           |                                                                                                                                                                                                                                                             |
|-----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 4.7       | Filtering: 1-6, 6 level                                                                                                                                                                                                                                     |
| 4.8       | Sound power: 2-6, 4 level                                                                                                                                                                                                                                   |
| 4.9       | Noise reduction: 0-4, 4 level                                                                                                                                                                                                                               |
| 4.10      | Smooth treatment: 0-4, 4 level                                                                                                                                                                                                                              |
| 4.11      | Frame correlation: 0-6, 6 level                                                                                                                                                                                                                             |
| 4.12      | Chromatography (Blood flow graph) : 0-37, 37 level                                                                                                                                                                                                          |
| 4.13      | Line density: Low-Middle-High, 3 level                                                                                                                                                                                                                      |
| 4.14      | Frequency: 4 level adjustable                                                                                                                                                                                                                               |
| 4.15      | Velocity: Minimum 0.4K, Maximum 40.5K<br>Convex probe: 0.4K-4.3K-38.5K<br>Linear probe: 0.4K-14.7K-39.0K<br>Trans-vaginal probe: 0.4K-7.8K-39.7K<br>Volume probe: 0.4K-4.2K-34.8K<br>Micro-convex probe: 0.4K-10.3K-40.5K<br>cardiac probe: 0.4K-7.8K-39.7K |
|           | PS: The frequency of the probe changes and the frequency value changes                                                                                                                                                                                      |
|           | PS: Frame rate changes with speed                                                                                                                                                                                                                           |
| <b>5:</b> | <b>Pulse wave Doppler (PW)</b>                                                                                                                                                                                                                              |
| 1         | Gain: 0-100, Step 2                                                                                                                                                                                                                                         |
| 5.2       | Spectrum envelope function: real time automatic spectrum envelope, manual spectrum envelope, and other modes.<br>The system automatically analyses and displays various data such as PSV, EDV, RI, PI, S/D, ACC, HR and so on.<br>Can wake up or close      |
| 5.3       | Sample volume: 0.5mm~30mm                                                                                                                                                                                                                                   |
| 5.4       | Blood angel: -75—75 degree, Step 5                                                                                                                                                                                                                          |
| 5.5       | False color: 0-67, 67 level                                                                                                                                                                                                                                 |
| 5.6       | Dynamic range: 20-40, 4 level                                                                                                                                                                                                                               |
| 5.7       | Filter: 0-9, 9 level                                                                                                                                                                                                                                        |
| 5.8       | Smooth treatment: 1-4, 4 level                                                                                                                                                                                                                              |
| 5.9       | Sound power: 2-5, 4 level                                                                                                                                                                                                                                   |
| 5.10      | Volume: 0-100, 10 level, Step 10                                                                                                                                                                                                                            |
| 5.11      | Audio filtering: 0-4, 4 level                                                                                                                                                                                                                               |
| 5.12      | Base line: -1.0~1.0,                                                                                                                                                                                                                                        |
| 5.13      | Grey map: 0-67, 67 level                                                                                                                                                                                                                                    |
| 5.14      | Scan velocity: 100-500, 6 level                                                                                                                                                                                                                             |
| 5.15      | PRF:Minimum 0.5K, Maximum 87.5K<br>Convex probe: 0.5K-4.3K-63.3K<br>Linear probe: 0.5K-14.5K-78.4K<br>Trans-vaginal probe: 0.5K-8.1K-78.4K                                                                                                                  |



DAWEI MEDICAL (JIANGSU) CORP., LTD.

|            |                                                                                                         |
|------------|---------------------------------------------------------------------------------------------------------|
|            | Volume probe: 0.5K-4.2K-53.8K<br>Micro-convex probe: 0.5K-10.3K-81.1K<br>cardiac probe: 0.5K-4.3K-87.5K |
| 5.16       | Frequency: 4 level                                                                                      |
|            | PS: The frequency of the probe changes and the PRF value changes                                        |
|            | PS: The frequency of the probe changes and the frequency value changes                                  |
| <b>6:</b>  | <b>Continuous Wave Doppler (CW)</b>                                                                     |
| 6.1        | Support probe: Cardiac probe                                                                            |
| 6.2        | Adjustment of B mode parameters is switchable                                                           |
| 6.3        | Gain: 0-100, Step 2                                                                                     |
| 6.4        | Sampling line position is adjustable                                                                    |
| 6.5        | PRF:0.9K~36.1K                                                                                          |
| 6.6        | Baseline: -1.0~1.0                                                                                      |
| 6.7        | Blood angel: -75~75 degree                                                                              |
| 6.8        | Grey map: 0-67                                                                                          |
| 6.9        | Scan velocity: 100-300                                                                                  |
| 6.10       | False color: 0-67                                                                                       |
| 6.11       | Dynamic range: 20-40                                                                                    |
| 6.12       | Filtering: 0-9, 9 level                                                                                 |
| 6.13       | Smooth treatment: 1-4                                                                                   |
| 6.14       | Frequency: 2.0MHz/2.3MHz/2.5MHz/3.0MHz, 4 level adjustable                                              |
| 6.15       | Sound power: 2-5                                                                                        |
| 6.16       | Volume: 0-100                                                                                           |
| 6.17       | Audio Filtering: 0-4                                                                                    |
| <b>☆7:</b> | <b>Anatomical M imaging</b>                                                                             |
| 7.1        | Support probe: Convex probe, Linear probe, Cardiac probe                                                |
| 7.2        | Adjustment of B mode parameters is switchable                                                           |
| 7.3        | Gain: 0-100, Step 2                                                                                     |
| 7.4        | M Sampling line angel is adjustable                                                                     |
| 7.5        | M Sampling line length is adjustable                                                                    |
| 7.6        | Sampling line: 3, Can be displayed or hidden separately                                                 |
| <b>☆8:</b> | <b>Color M mode (MC)</b>                                                                                |
| 8.1        | Adjustment of B mode parameters is switchable                                                           |
| 8.2        | Gain: 0-100, Step2                                                                                      |
| 8.3        | MC Sampling line angel is adjustable                                                                    |



DAWEI MEDICAL (JIANGSU) CORP., LTD.

|             |                                                                  |
|-------------|------------------------------------------------------------------|
| 8.4         | MC Sampling line length is adjustable                            |
| 8.5         | Frequency: 4 level                                               |
| 8.6         | Sampling number: 6-24                                            |
| 8.7         | Speed through: 0-8, 8 level                                      |
| 8.8         | Scan velocity: 150-500                                           |
| 8.9         | Frame correlation: 0-6, 6 level                                  |
| 8.10        | Filtering: 1-6, 6 level                                          |
| 8.11        | Blood flow preferred: 0-8, 8 level                               |
| 8.12        | Smooth treatment: 0-4, 4 level                                   |
| 8.13        | Map: 0-37, 37 level                                              |
| <b>☆9:</b>  | <b>Elastography</b>                                              |
| 9.1         | Adjustment of B mode parameters is switchable                    |
| 9.2         | Gain: 0-100, Step 2                                              |
| 9.3         | B/E, Double real-time display on the same screen                 |
| 9.4         | Probe displacement curve display: Up/Down                        |
| 9.5         | Pressure indicator bar display                                   |
| 9.6         | Frequency: 8-9 level, Adjustable; According to the probe display |
| 9.7         | Noise reduction: 0-2, 2 level                                    |
| 9.8         | Frame correlation: 0-3, 3 level                                  |
| 9.9         | Comparison: 0-13, 13 level                                       |
| 9.10        | False color: 0-3, 3 level                                        |
| 9.11        | Don't support cardiac probe                                      |
| <b>☆10:</b> | <b>Tissue Doppler imaging (TDI)</b>                              |
| 10.1        | Support probe: Cardiac probe                                     |
| 10.2        | Adjustment of B mode parameters is switchable                    |
| 10.3        | Gain: 0-100, step 2                                              |
| 10.4        | ROI area adjustable                                              |
| 10.5        | Sampling number: 6-24                                            |
| 10.6        | Velocity: 0.4K-8.0K                                              |
| 10.7        | Frame correlation: 0-6, 6 level                                  |
| 10.8        | Tissue preferred: 0-7, 7 level                                   |
| 10.9        | Frequency: 2.0MHz/2.3MHz/2.5MHz/3.0MHz                           |
| 10.10       | Support color reversal                                           |
| <b>☆11:</b> | <b>Strain rate imaging</b>                                       |
| 11.1        | Support probe: Cardiac probe                                     |



DAWEI MEDICAL (JIANGSU) CORP., LTD.

|             |                                               |
|-------------|-----------------------------------------------|
| 11.2        | Adjustment of B mode parameters is switchable |
| 11.3        | ROI area adjustable                           |
| 11.4        | Gain: 0-100, Step 2                           |
| 11.5        | Sampling number: 6-24, 6 level                |
| 11.6        | Axial average: 1-4, 4 level                   |
| 11.7        | Velocity: 0.4K-8K                             |
| 11.8        | Frame correlation: 0-6, 6 level               |
| 11.9        | Tissue optimization: 0-7, 7 level             |
| <b>☆12:</b> | <b>Panoramic imaging</b>                      |
| 12.1        | Support probe: Linear probe                   |
| 12.2        | Speckle Reduction: 0-5, 5 level               |
| <b>☆13:</b> | <b>Deflection imaging</b>                     |
| 13.1        | Support probe: Linear probe                   |
| 13.2        | Adjustment of B mode parameters is switchable |
| 13.3        | Deflection angel: 8 level                     |
| 13.4        | Speckle reduction: 0-5, 5 level               |
| 13.5        | Dynamic rate: 30-180, Step 5                  |
| 13.6        | Line density: low-middle-high, 3 level        |
| 13.7        | Frame Correlation: 0-4, 4 level               |
| 13.8        | False color: 0-67, 67 level                   |
| 13.9        | Image style: Soft-Comparison, 2 level         |
| 13.10       | Noise reduction: 0-5, 5 level                 |
| 13.11       | Edge Enhancement: 0-5, 5 level                |
| 13.12       | Sound power: 2-10, 8 level                    |
| 13.13       | Grey map: 0-67, 67 level                      |
| <b>☆14:</b> | <b>Trapezoidal imaging</b>                    |
| 14.1        | Probe support: linear probe                   |
| 14.2        | Adjustment of B mode parameters is switchable |
| 14.3        | Deflection angel: 8 level                     |
| 14.4        | Speckle reduction: 0-5, 5 level               |
| 14.5        | Dynamic rate: 30-180, Step 5                  |
| 14.6        | Line density: low-middle-high, 3 level        |
| 14.7        | Frame Correlation: 0-4, 4 level               |
| 14.8        | False color: 0-67, 67 level                   |
| 14.9        | Image style: Soft-Comparison, 2 level         |

|            |                                                                                   |
|------------|-----------------------------------------------------------------------------------|
| 14.10      | Noise reduction: 0-5, 5 level                                                     |
| 14.11      | Edge Enhancement: 0-5, 5 level                                                    |
| 14.12      | Sound power: 2-10, 8 level                                                        |
| 14.13      | Grey map: 0-67, 67 level                                                          |
| 14.14      | Space Synthesis: 0-2, 2 level                                                     |
| <b>☆15</b> | <b>Freehand 3D imaging</b>                                                        |
| 15.1       | Support probe: convex probe, linear probe                                         |
| 15.2       | Display model: 4 pictures                                                         |
| 15.3       | Image Rotation X/Y/Z Axis                                                         |
| 15.4       | Multi-slice Visibility                                                            |
| <b>16</b>  | <b>Real-time 4D imaging</b>                                                       |
| 16.1       | Support probe: 4D volume probe                                                    |
| 16.2       | Adjustment of B mode parameters is switchable                                     |
| 16.3       | Gain: 0-100, Step 2                                                               |
| 16.4       | Display model: one image、two images、four images                                   |
| 16.5       | Image Rotation: X/Y/Z Axis                                                        |
| 16.6       | Multi-slice Visibility                                                            |
| 16.7       | Light&Shade inversion                                                             |
| 16.8       | Smooth: 0-4, 4 level                                                              |
| 16.9       | Threshold level: 0-129, Step 3                                                    |
| 16.10      | Transparency: 1-509, Step 10                                                      |
| 16.11      | Render type: 4 kinds, Surface、maximum、minimum、perspective                         |
| <b>17:</b> | <b>Extended Imaging</b>                                                           |
| 17.1       | Gain: 0—100, Step 2                                                               |
| 17.2       | TGC: 8 segment adjustable                                                         |
| 17.3       | Maximum focus point: $\geq 7$ , which can be moved throughout the whole process.  |
| 17.4       | Speckle reduction: 0-5, 5 level                                                   |
| 17.5       | Space Compound: 0-2, 2 level (Linear probe: 3 level, don't support cardiac probe) |
| 17.6       | Dynamic range: 30-180, 35 level, Step 5                                           |
| 17.7       | Line density: Low、Middle、High, 3 level                                            |
| 17.8       | Frame correlation: 0-4, 4 level                                                   |
| 17.9       | Noise reduction: 0-5, 5 level                                                     |
| 17.10      | Edge enhancement: 0-5, 5 level                                                    |
| 17.11      | Sound power: 2-10, 9 level                                                        |
| 17.12      | Grey map: 0-67, 67 level                                                          |

|            |                                                                                                                                                                                                  |
|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 17.13      | False color: 0-67, 67 level                                                                                                                                                                      |
| 17.14      | Image style: Soft-Comparison, 2 level                                                                                                                                                            |
| 17.15      | Extended level: Maximum 72 level<br>Convex probe: 9 level<br>Trans-vaginal probe: 72 level<br>Micro-convex probe: 29 level<br>Cardiac probe: 40 level<br>4D Volume probe: 17 level               |
|            | PS: The screen has real-time display of voice power, probe frequency, dynamic range, pseudo color, gray scale and other 11 parameters can be adjusted                                            |
|            | PS: When the probe scan range reaches the maximum, the space synthesized is 0.                                                                                                                   |
| <b>18:</b> | <b>Measurement and analysis function:</b>                                                                                                                                                        |
| 18.1       | General measurement: Distance, area, ellipse, cross line, angle, distance ratio, volume, Volume (ellipse), area ratio, diameter, joint angle                                                     |
| 18.2       | Cardiac: Automatic spectrum envelope、LV、Main Pulmonary artery diameter、RVEDd、RVEDs、LVM、LAV、HR、MVF、AO、AR、LVOT、TVF、Pulmonic valve、Pulmonary vein、RV、Doppler fetal heart sound、LVET、LVM、LVMI、AV     |
| 18.3       | Vascular: carotid intima (IMT), length stenosis ratio, area stenosis ratio, IMT (back wall), IMT (front wall)                                                                                    |
| 18.4       | OB: Fetal routine、AFI、TW、GS、CRL、OFD、HL、ulna、NT、Fibula、Nbone1、Radial、Tibia                                                                                                                        |
| 18.5       | GYN: uterus、cervix、corpus uteri/cervix uterus、left ovarian vein、right ovarian vein、dominant follicle、intima thickness                                                                            |
| 18.6       | Urology: prostate、residual urine、left kidney、right kidney、left suprarenal vein、right suprarenal、left testis、right testis、left seminal vesicle、right seminal vesicle                              |
| 18.7       | Abdomen: liver、CHD、portal vein diameter、cholecyst、CBD、pancreas、spleen、Internal diameter of abdominal aorta、kidney                                                                                |
| 18.8       | Small parts: Thyroid                                                                                                                                                                             |
| 18.9       | Software package: Measurement package、Software package、Medical records management software package                                                                                               |
| <b>19:</b> | <b>Graphic and text management system</b>                                                                                                                                                        |
| 19.1       | Host build in hard disk (SSD 128GB), Start fast and stable, RAM: 8GB                                                                                                                             |
| 19.2       | Movie playback: $\geq 1200$ frames                                                                                                                                                               |
| 19.3       | Internal file information management system: can record patient number, name, check number, check date and so on, and can be searched and managed by numbering, checking number, name and so on. |
| 19.4       | Type of report is 16                                                                                                                                                                             |
| 19.5       | One key fast report graphic and text management                                                                                                                                                  |
| <b>20:</b> | <b>Interface</b>                                                                                                                                                                                 |
| 20.1       | USB interface: 2                                                                                                                                                                                 |

|            |                                                                                                                                                     |
|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| 20.2       | HDMI interface: 1                                                                                                                                   |
| 20.3       | RJ-45 interface: 1                                                                                                                                  |
| 20.4       | Grounding wire interface: 1                                                                                                                         |
| <b>21:</b> | <b>Configuration</b>                                                                                                                                |
| 21.1       | Portable type full digital color Doppler ultrasound diagnostic system                                                                               |
| 21.2       | Probe: convex array probe (standard), linear probe (optional), Trans-vaginal probe (optional), cardiac probe (optional), 4D volume probe (optional) |
| 21.3       | ≥13 quick adjusting knobs                                                                                                                           |
| <b>22:</b> | <b>Technology, after-sales service and other requirements</b>                                                                                       |
| 22.1       | After acceptance, the warranty is free for two years (Provide manufacturer warranty certificate)                                                    |
| 22.2       | Manufacturer has ISO13485 certification and EU CE certification.                                                                                    |
|            |                                                                                                                                                     |

## Image Gallery



