

# HD10 Vet

## Veterinary Patient Monitor



### Short boot time

10s to enter working state



### Reliable

Extremely strong anti-interference ability, stable performance



### HD display

12.1" wide screen, high definition display



### Quick response

Millisecond touch response



ECG Full Screen



Large Font



Standard Interface



No-fan Design



Night Mode



3 Levels Adjustable Volume



Shortcut Menu



Thermal Recorder

### ECG:

Advanced Ag/AgCl electroshock measurement method to provide physicians with accurate ECG measurements. More than six kinds of gain options, full-screen mode without overlapping each other.

### NIBP:

Anti-interference and antiweak signal algorithm designed for the characteristics of animal blood flow, enables fast measurement of accurate blood pressure values.

### SpO2:

Newly developed animal blood oxygen module with super resistance to motion and weak perfusion capability, applied to different size animals and quickly measure accurate values.

### CO2:

A high sampling rate CO2 module that provides accurate and clear CO2 tracings at 160 breaths/min. The ultra-low flow rate design allows for accurate sampling at 50 ml/min, ideal for all types of diseased animals.



Nurse Call



USB Output



Defib Proof & Sync



HDMI Video Output  
(Optional)



Alarm Light Strip



2500  
NIBP Measurement  
Review



30 days  
Review Storage



300min  
Continuous Working Time



7000  
Event Log




10s  
Boot Time



Strict Quality Control



Plug-in Accessories



Touch Screen  
(Optional)